

Minimally Invasive Surgery

Surgical Skills



Lesson Objectives:

1. Describe the preparation of the patient for minimally invasive surgery (MIS)
2. Describe the function of each component of the imaging equipment used in MIS
3. Describe the surgical technique used for insufflation in laparoscopy
4. Describe the trocar-cannula system used in MIS
5. Describe the specific electrosurgical risks of direct and capacitive coupling
6. Describe the structure and function of a flexible endoscope

Minimally Invasive Surgery (MIS)

Surgical Techniques performed through small incisions or natural body openings.


- **Goals:**

- Minimize trauma to the body.
- Reduce postoperative pain.
- Shorten recovery time.
- Lower risk of complications compared to open surgery.

- **Examples:**

- Laparoscopic surgery.
- Arthroscopy.
- Endoscopic surgery.
- Robotic-assisted surgery.

Flow of Ideas




Technology and Techniques in Minimally Invasive Surgery

The imaging system must produce a clear view of the anatomy in real time to achieve the surgical goals while supporting safe surgery.

Imaging System includes:

- Surgical telescope (or endoscope)
- Camera head
- Camera control unit (CCU)
- Light source
- Light cable
- Video cables
- Monitor screen
- Image management system



Minimally Invasive Surgery Instruments

- **Design Features:**
 - Handles and fulcrums at a distance from working end.
 - Small hinges, springs, and valves.
 - Rotational design for increased maneuverability.
- **Common Instruments:**
 - Retractors: Extend from shaft tip, flare out or curve for tissue manipulation.
 - Scissors: Straight, curved, and hooked configurations.
 - Grasping Instruments: Clamps and forceps, with atraumatic or penetrating tips.



Energy Sources in Minimally Invasive Surgery

Electrosurgery

Ultrasonic Devices


Lasers

Components of Electrosurgical Unit (ESU) Circuit

- ESU
 - Generates power
- Active electrode
 - Hand piece/pencil
- Patient
- Dispersive electrode (grounding pad)

Flow of Electricity:

- Generator
- Active electrode
- Surgical site
 - Heat (Hemostasis)
- Patient
- Dispersive electrode (G- pad)
 - Generator



Electrosurgical Risks

- Insulation Failure: Risk of burns due to damaged insulation.
- Direct and Capacitative Coupling: Burns from contact with conductive instruments.
- **Risk Reduction:**
 - Active electrode monitoring (AEM).
 - Regular inspection for damage.
 - Avoidance of hybrid cannulas.
 - Monitoring power settings.
 - Proper reprocessing of instruments.
 - Preventing cord damage and wet surface contact.

Ultrasonic Energy

- Ultrasonic technology in MIS for coagulation and cutting.
- Coagulates tissue by creating cool coagulum at cellular level.
- Minimal heat and no electrical energy involved.
- Systems include SonoSurg (Olympus America) and Harmonic scalpel (Ethicon, Johnson & Johnson).

Lasers

- Utilized in minimally invasive and endoscopic surgery.
- Lasing fiber introduced through cannula or telescope.
- Advances in digital technology enhance precision.
- Hazards include limited space and blind areas.

Four Interactions of Laser and Tissue

- **Absorption**
 - Tissue destruction
 - Heat, pressure, rupture
- **Transmission**
 - Through media
 - Water, aqueous humor



- **Reflection**
 - Mirror
 - Accidental
 - Inadvertent burn
 - Non-reflective Inst
- **Scattering**
 - Broken up & dispersed

Patient Preparation for MIS

Patient Positioning

- Determined by Surgical Site: Varied positions for different procedures.
- Attention to Safety:
 - Immediate adjustment for intraoperative changes.
 - Prevention of patient sliding with appropriate devices.

Skin Prep and Draping

- Preparation for Conversion: Skin prep as for open surgery.
- Extended Draping: Matched to potential open incision area.

Surgical Setup for MIS

Organized Instrument Placement

Separate MIS instruments to prevent damage by heavier equipment.

Securement of Cords and Tubing

Use nonpenetrating clamps to secure cords and tubing to drapes.

Passing Instruments

Orient instruments with working tip downward when passing to surgeon.

Lighting Considerations

Dim or power off room lights during MIS

Converting to an Open Case

- Any minimally invasive surgery has the potential to become open case
- Perform this rapidly
- Communicate with circulator
- Quickly exchange instruments
- Count added instruments and sponges
- Initial counts should include all items that would need to be counted in the open procedure

Insufflation Used in Laparoscopy

- **Insufflation Overview**

- Process: Fills abdominal/thoracic cavities with CO₂ for visibility and instrument entry.
- Benefits: CO₂ is nontoxic, absorbable, nonflammable.
- Warming: Maintains core temp, prevents lens fogging.

- **Creating Pneumoperitoneum**

- Method: Surgeon inflates abdomen with CO₂ via Veress needle/trocar.
- Veress Needle: Springs retract upon resistance.
- Control: Adjust for steady gas delivery.
- Trocar Insertion: After abdomen inflation.

Watch the "Basic Laparoscopy: Access and Port placement" Video for an overview

Basic Laparoscopy: Access and Port placement Video



Basic Laparoscopy: Access and Port placement Video

- Summary of Video:
 - Access with insufflation needle – fill abdomen with CO₂ (pneumoperitoneum)
 - Access may also be with an optical trocar port

Risks with MIS

- **Insufflation failure**

- Can happen if pressure is lost in body cavity. Can be dangerous based on instruments in the body cavity at the time
- Insufflation can also put pressure on surrounding body parts – this is notable in thoracoscopy, as this pressure may affect respiratory and heart function. Pressure and flow should be balanced with the patient's ability to tolerate

- Insulation Failure – Many electrocautery MIS instruments have a protective sheathing
- Direct coupling and capacitive coupling

Risk Reduction

- Active electrode monitoring
- Inspect instruments
- Correct cannula use
- Power setting problems
- Reprocessing
- Cord placement
- Prevent kinks and knots
- Avoid wet surfaces with cords

**Watch the "Laparoscopic Instruments" Video for an
overview**

Laparoscopic Instruments Video

[Click Here](#) to watch the video!

Laparoscopic Instruments Video

Summary of Video:

- Parts of instrument: Tip, Shaft, Handle
- Instrument is designed based on use
- Instrument works in a mirrored function – the port is a fulcrum
- Triangulation: Port placement and working area is triangulated

Principles of Flexible Endoscopy

- Viewing inside body passages or organs
- Introduced through a natural opening
- Advanced carefully
- Used for Biopsies and examination



Watch the "What is Endoscopy" Video

What is Endoscopy Video



What is Endoscopy Video

- Summary of Video:
 - Flexible scope/camera inserted through a natural body opening (Mouth, Anus, etc)
 - Can be done with sedation in outpatient setting

Types of Endoscopes

Flexible Endoscopes

- Structure: Head and insertion tube.
- Components: Control head, fiber-optic light channel, instrument channel.
- Functionality: Illumination, flexing mechanism, rotational views.

Echoendoscopes

- Structure: Similar to standard endoscope.
- Additional Feature: Ultrasound transducer at the tip.
- Functionality: Transmitting and receiving ultrasound signals for imaging.

Types of Endoscopes




Capsule Endoscope

- Structure: Small capsule device.
- Components: Camera, lens, LED, battery.
- Functionality: Swallowed by patient, transmits images to sensors.

Imaging System

- Similarity to MIS: Camera control unit and documentation system.
- Usage in MIS: Assisted procedures for enhanced visibility.
- Combined Procedures: Abdomen and gastrointestinal tract enhanced with both technologies.

Endoscopes



Read Chapter 20 from the E-Book

Read **Chapter 20** from your E-Book to pass the upcoming quiz from **Surgical Technology - Elsevier eBook on VitalSource, 8th Edition**.

[Click Here](#) to access Chapter 20!

Thank you!

Get ready for your quiz and rest of the activities now. Best of luck!



Congratulations!

Lesson 20 is complete.